

Product datasheet for TA600304

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CBX8 Mouse Monoclonal Capture Antibody [Clone ID: OTI2C8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C8
Applications: LMNX

Recommended Dilution: 1:100 - 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CBX8 (NM_020649) produced in HEK293T

cell

Formulation: Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

Concentration: 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

Stability: Stable for 12 months from date of receipt.

Gene Name: chromobox 8

Database Link: NP 065700

Entrez Gene 30951 MouseEntrez Gene 303731 RatEntrez Gene 57332 Human

Q9HC52

Synonyms: PC3; RC1


Matched ELISA Pair: TA700302

Protein Families: Transcription Factors






Product images:



CBX8 Luminex Elisa with 2C8 Capture (TA600304) and 2A5 Detection ([TA700302]) Antibodies. Substrate used: Recombinant Human CBX8 ([LY412158])



CBX8 Luminex Elisa with 2C8 Capture (TA600304) and 3D1 Detection ([TA700303]) Antibodies. Substrate used: Recombinant Human CBX8 ([LY412158])