

Product datasheet for **TA600313**

CD1C Mouse Monoclonal Capture Antibody [Clone ID: OTI1F4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1F4
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD1C (NM_001765) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD1c molecule
Database Link:	NP_001756 Entrez Gene 911 Human P29017
Synonyms:	BDCA1; CD1; CD1A; R7
Matched ELISA Pair:	TA700313, TA700314
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage




[View online »](#)


Product images:




CD1C Elisa with 1F4 Capture (TA600313) and 2F4 Detection ([TA700313]) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Elisa with 1F4 Capture (TA600313) and 5A8 Detection ([TA700314]) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 1F4 Capture (TA600313) and 2F4 Detection ([TA700313]) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 1F4 Capture (TA600313) and 5A8 Detection ([TA700314]) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])