

## **Product datasheet for TA600334**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## MEK1 (MAP2K1) Mouse Monoclonal Capture Antibody [Clone ID: OTI3E6]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3E6
Applications: ELISA

Recommended Dilution: ELISA 1:100-1:1000

**Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MAP2K1 (NM\_002755) produced in

HEK293T cell.

**Formulation:** Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

Concentration: 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

Stability: Stable for 12 months from date of receipt.

Gene Name: mitogen-activated protein kinase kinase 1

Database Link: NP 002746

Entrez Gene 26395 MouseEntrez Gene 170851 RatEntrez Gene 5604 Human

Q02750

Synonyms: CFC3; MAPKK1; MEK1; MKK1; PRKMK1

Matched ELISA Pair: TA700334

**Protein Families:** Druggable Genome, Protein Kinase







### **Protein Pathways:**

Acute myeloid leukemia, B cell receptor signaling pathway, Bladder cancer, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Dorso-ventral axis formation, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Gap junction, Glioma, GnRH signaling pathway, Insulin signaling pathway, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Melanoma, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Oocyte meiosis, Pancreatic cancer, Pathways in cancer, Prion diseases, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, T cell receptor signaling pathway, Thyroid cancer, Toll-like receptor signaling pathway, Vascular smooth muscle contraction, VEGF signaling pathway

# **Product images:**



MAP2K1 Elisa with 3E6 Capture (TA600334) and 7B1 Detection ([TA700334]) Antibodies. Substrate used: Recombinant Human MAP2K1 ([TP318460])



MAP2K1 Luminex Elisa with 3E6 Capture (TA600334) and 7B1 Detection ([TA700334]) Antibodies. Substrate used: Recombinant Human MAP2K1 ([TP318460])