

Product datasheet for **TA600367**


FGFR2 Mouse Monoclonal Capture Antibody [Clone ID: OTI5F7]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI5F7
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FGFR2 (NM_000141) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fibroblast growth factor receptor 2
Database Link:	NP_000132 Entrez Gene 14183 Mouse Entrez Gene 25022 Rat Entrez Gene 2263 Human P21802
Synonyms:	BBDS; BEK; BFR-1; CD332; CEK3; CFD1; ECT1; JWS; K-SAM; KGFR; TK14; TK25
Matched ELISA Pair:	TA700365
Protein Families:	Druggable Genome, Protein Kinase, Secreted Protein, Transmembrane
Protein Pathways:	Endocytosis, MAPK signaling pathway, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton



[View online »](#)

Product images:


FGFR2 ELISA with 5F7 Capture (TA600367) and 3B3 Detection ([TA700365]) Antibodies. Substrate used: Recombinant Human FGFR2 ([TP317098])



FGFR2 Luminex Elisa with 5F7 Capture (TA600367) and 3B3 Detection ([TA700365]) Antibodies. Substrate used: Recombinant Human FGFR2 ([TP317098])