

## Product datasheet for **TA600369**

### FGFR2 Mouse Monoclonal Capture Antibody [Clone ID: OTI3F8]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3F8
Applications:	LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FGFR2 (NM_000141) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	fibroblast growth factor receptor 2
Database Link:	<a href="#">NP_000132</a> <a href="#">Entrez Gene 14183 Mouse</a> <a href="#">Entrez Gene 25022 Rat</a> <a href="#">Entrez Gene 2263 Human</a> <a href="#">P21802</a>
Synonyms:	BBDS; BEK; BFR-1; CD332; CEK3; CFD1; ECT1; JWS; K-SAM; KGFR; TK14; TK25
Matched ELISA Pair:	TA700366
Protein Families:	Druggable Genome, Protein Kinase, Secreted Protein, Transmembrane
Protein Pathways:	Endocytosis, MAPK signaling pathway, Pathways in cancer, Prostate cancer, Regulation of actin cytoskeleton



[View online »](#)

## Product images:



FGFR2 Luminex Elisa with 3F8 Capture (TA600369) and 1H5 Detection ([TA700366]) Antibodies. Substrate used: Recombinant Human FGFR2 ([TP317098])