

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA600373

CD33 Mouse Monoclonal Capture Antibody [Clone ID: OTI1E4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1E4
Applications:	ELISA, LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD33 (NM_001772) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD33 molecule
Database Link:	<u>NP_001763</u> <u>Entrez Gene 945 Human</u> <u>P20138</u>
Synonyms:	p67; SIGLEC-3; SIGLEC3
Matched ELISA Pair:	ТА700372, ТА700373
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


Product images:



CD33 ELISA with 1E4 Capture (TA600373) and 1B10 Detection ([TA700372]) Antibodies. Substrate used: Recombinant Human CD33 ([TP307023])

CD33 ELISA with 1E4 Capture (TA600373) and 2F8 Detection ([TA700373]) Antibodies. Substrate used: Recombinant Human CD33 ([TP307023])

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



CD33 Luminex Elisa with 1E4 Capture (TA600373) and 1B10 Detection ([TA700372]) Antibodies. Substrate used: Recombinant Human CD33 ([TP307023])

CD33 Luminex Elisa with 1E4 Capture (TA600373) and 2F8 Detection ([TA700373]) Antibodies. Substrate used: Recombinant Human CD33 ([TP307023])

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US