

## Product datasheet for **TA600382**

### **p63 (TP63) Mouse Monoclonal Capture Antibody [Clone ID: OTI5G7]**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI5G7
<b>Applications:</b>	ELISA, LMNX
<b>Recommended Dilution:</b>	1:100 - 1:1000
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human TP63 (NM_003722) produced in HEK293T cell.
<b>Formulation:</b>	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
<b>Concentration:</b>	0.5 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store at -20°C. Avoid repeat freeze/thaw cycles.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	tumor protein p63
<b>Database Link:</b>	<a href="#">NP_003713</a> <a href="#">Entrez Gene 246334 Rat</a> <a href="#">Entrez Gene 8626 Human</a> <a href="#">Q9H3D4</a>
<b>Synonyms:</b>	AIS; B(p51A); B(p51B); EEC3; KET; LMS; NBP; OFC8; p40; p51; p53CP; p63; p73H; p73L; RHS; SHFM4
<b>Matched ELISA Pair:</b>	TA700382
<b>Protein Families:</b>	Druggable Genome, Transcription Factors




[View online »](#)

Product images:



TP63 ELISA with 5G7 Capture (TA600382) and 5F3 Detection ([TA700382]) Antibodies. Substrate used: Recombinant Human TP63 ([TP308013])



TP63 Luminex Elisa with 5G7 Capture (TA600382) and 5F3 Detection ([TA700382]) Antibodies. Substrate used: Recombinant Human TP63 ([TP308013])