

Product datasheet for **TA600434**


Frataxin (FXN) Mouse Monoclonal Capture Antibody [Clone ID: OTI4G4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4G4
Applications:	LMNX
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human FXN (NM_000144) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	frataxin
Database Link:	NP_000135 Entrez Gene 2395 Human Q16595
Synonyms:	CyaY; FA; FARR; FRDA; X25
Matched ELISA Pair:	TA700050
Protein Families:	Druggable Genome



[View online »](#)

Product images:

FXN Luminex ELISA with 4G4 Capture (TA600434) and Anti-Tag 1B12 Detection ([TA700050]) Antibodies. Substrate used: Recombinant Human FXN ([TP304880])