

Product datasheet for **TA600466**

CD5 Mouse Monoclonal Capture Antibody [Clone ID: OTI5D4]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5D4
Applications:	ELISA
Recommended Dilution:	1:100 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD5 (NM_014207) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD5 molecule
Database Link:	NP_055022 Entrez Gene 921 Human P06127
Synonyms:	LEU1; T1
Matched ELISA Pair:	TA700466
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage



[View online »](#)

Product images:



CD5 ELISA with 5D4 Capture (TA600466) and 10H3 Detection ([TA700466]) Antibodies. Substrate used: Recombinant Human CD5 ([TP306494])