

Product datasheet for TA600491

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

B7-1 (CD80) Mouse Monoclonal Capture Antibody [Clone ID: OTI2E5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2E5
Applications: LMNX

Recommended Dilution: 1:100 - 1:1000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CD80 (NM_005191) produced in HEK293T

cell

Formulation: Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose

Concentration: 0.5 mg/ml

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

Stability: Stable for 12 months from date of receipt.

Gene Name: CD80 molecule

Database Link: NP 005182

Entrez Gene 941 Human

P33681

Synonyms: B7; B7-1; B7.1; BB1; CD28LG; CD28LG1; LAB7

Matched ELISA Pair: TA700489

Protein Families: Druggable Genome, Transcription Factors, Transmembrane

Protein Pathways: Allograft rejection, Autoimmune thyroid disease, Cell adhesion molecules (CAMs), Graft-


versus-host disease, Systemic lupus erythematosus, Toll-like receptor signaling pathway,

Type I diabetes mellitus, Viral myocarditis





Product images:



CD80 Luminex ELISA with 2E5 Capture (TA600491) and 2B11 Detection ([TA700489]) Antibodies. Substrate used: Recombinant Human CD80 ([TP306540])