

Product datasheet for **TA600525**

IKZF3 Mouse Monoclonal Capture Antibody [Clone ID: OTI2B9]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2B9
Applications:	ELISA
Recommended Dilution:	ELISA 1:100-1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human IKZF3 (NM_012481) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) containing 0.05% sodium azide and up to 5% trehalose
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	IKAROS family zinc finger 3
Database Link:	NP_036613 Entrez Gene 22780 Mouse Entrez Gene 303511 Rat Entrez Gene 22806 Human Q9UKT9
Synonyms:	AIO; AIOLOS; ZNFN1A3
Matched ELISA Pair:	TA700525




[View online »](#)

Product images:



IKZF3 ELISA with 2B9 Capture (TA600525) and 5H2 Detection ([TA700525]) Antibodies. Substrate used: Recombinant Human IKZF3 ([TP307547])



IKZF3 Luminex with 2B9 Capture (TA600525) and 5H2 Detection ([TA700525]) Antibodies. Substrate used: Recombinant Human IKZF3 ([TP307547])