

Product datasheet for TA600555

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

IL17 (IL17A) Mouse Monoclonal Capture Antibody [Clone ID: OTI3D5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3D5

Applications: ELISA, LMNX

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human IL-17A (NP_002181) produced in 293 cell.

Formulation: PBS (pH7.3) containing 0.05% sodium azide

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated Storage: Store at -20°C.

Avoid repeat freeze/thaw cycles.

Predicted Protein Size: 17.5 kDa

Gene Name: interleukin 17A

Database Link: NP 002181

Entrez Gene 3605 Human

Q16552

Background: The protein encoded by this gene is a proinflammatory cytokine produced by activated T

cells. This cytokine regulates the activities of NF-kappaB and mitogen-activated protein kinases. This cytokine can stimulate the expression of IL6 and cyclooxygenase-2 (PTGS2/COX-2), as well as enhance the production of nitric oxide (NO). High levels of this cytokine are associated with several chronic inflammatory diseases including rheumatoid arthritis,

psoriasis and multiple sclerosis. [provided by RefSeq, Jul 2008]

Synonyms: CTLA-8; CTLA8; IL-17; IL-17A; IL17; ILA17

Matched ELISA Pair: TA600556, TA700556







Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Cytokine-cytokine receptor interaction

Product images:



IL17A ELISA with 3D5 Capture (TA600555) and 1E4 Detection ([TA700556])) Antibodies. Substrate used: Recombinant Human IL17A (KEK293 cells expressed protein)



IL17A Luminex ELISA with 3D5 Capture (TA600555) and 1E4 Detection ([TA700556])) Antibodies. Substrate used: Recombinant Human IL17A (KEK293 cells expressed protein)