

## Product datasheet for **TA600561**

### **gag Rabbit Monoclonal Capture Antibody [Clone ID: OTIR5D12]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTIR5D12
Applications:	ELISA
Reactivity:	HIV
Host:	Rabbit
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Recombinant protein fragment corresponding to amino acids 1-231 of HIV P24 produced in E.coli.
Formulation:	PBS (PH 7.3) containing 0.05% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C. Avoid repeat freeze/thaw cycles.
Predicted Protein Size:	25.58kDa
Matched ELISA Pair:	TA700561



[View online »](#)

Product images:



HIV P24 ELISA with 5D12 Capture (TA600561) and 2C5 Detection ([TA700561]) Antibodies. Substrate used: Recombinant HIV P24 (E.coli expressed protein)