

Product datasheet for **TA700032**

BIRC5 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI3G3]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3G3
Applications:	ELISA, LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human BIRC5 (NP_001159) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	baculoviral IAP repeat containing 5
Database Link:	NP_001159 Entrez Gene 332 Human O15392
Synonyms:	API4; EPR-1
Matched ELISA Pair:	TA600032
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Colorectal cancer, Pathways in cancer




[View online »](#)

Product images:



BIRC5 ELISA with 2D9 Capture ([TA600032]) and 3G3 Detection (TA700032) Antibodies. Substrate used: Recombinant Human BIRC5 ([TP305935])



BIRC5 Luminex ELISA with 2D9 Capture ([TA600032]) and 3G3 Detection (TA700032) Antibodies. Substrate used: Recombinant Human BIRC5 ([TP305935])