

Product datasheet for TA700040

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GFAP Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI2F7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F7

Applications: ELISA, LMNX **Recommended Dilution:** 1:250 - 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human GFAP (NM_002055) produced in HEK293T

cell.

Formulation: Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol

Concentration: 0.5 mg/ml **Conjugation:** Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: glial fibrillary acidic protein

Database Link: NP 002046

Entrez Gene 14580 MouseEntrez Gene 24387 RatEntrez Gene 2670 Human

P14136

Synonyms: ALXDRD


Matched ELISA Pair: TA500337, TA600040

Protein Families: ES Cell Differentiation/IPS






Product images:



GFAP ELISA with 4G8 Capture ([TA500337]) and 2F7 Detection (TA700040) Antibodies. Substrate used: Recombinant Human GFAP ([TP304548])



GFAP Luminex with 5D7 Capture ([TA600040]) and 2F7 Detection (TA700040) Antibodies. Substrate used: Recombinant Human GFAP ([TP304548])