

Product datasheet for **TA700050**

DDK Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: 1B12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1B12
Applications:	ELISA
Recommended Dilution:	ELISA 1:250-1:1000
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Matched ELISA Pair:	TA600098, TA600144, TA600251, TA600252, TA600253, TA600254, TA600255, TA600256, TA600257, TA600258, TA600259, TA600260, TA600309, TA600310, TA600324, TA600325, TA600406, TA600407, TA600416, TA600417, TA600432, TA600433, TA600434, TA600486, TA600487, TA600488, TA600498, TA600499, TA600514, TA600528, TA600535, TA600536, TA600537



[View online »](#)