

Product datasheet for **TA700125**

SOX17 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI3H5]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3H5
Applications:	LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 177-414 of human SOX17 (NP_071899) produced in E.coli.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SRY-box transcription factor 17
Database Link:	NP_071899 Entrez Gene 20671 Mouse Entrez Gene 312936 Rat Entrez Gene 64321 Human Q9H6I2
Synonyms:	VUR3
Matched ELISA Pair:	TA600124
Protein Families:	Transcription Factors
Protein Pathways:	Wnt signaling pathway




[View online »](#)

Product images:



SOX17 Luminex with 2G8 Capture ([TA600124]) and 3H5 Detection (TA700125) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 Luminex with 5G6 Capture ([TA600126]) and 3H5 Detection (TA700125) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])