

## Product datasheet for **TA700126**

### **SOX17 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI5G6]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI5G6
Applications:	ELISA, LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant protein expressed in E.coli corresponding to amino acids 177-414 of human SOX17
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	SRY-box transcription factor 17
Database Link:	<a href="#">NP_071899</a> <a href="#">Entrez Gene 20671 MouseEntrez Gene 312936 RatEntrez Gene 64321 Human Q9H6I2</a>
Synonyms:	VUR3
Matched ELISA Pair:	TA600123
Protein Families:	Transcription Factors
Protein Pathways:	Wnt signaling pathway




[View online »](#)


**Product images:**




SOX17 ELISA with 2F9 Capture ([TA600123]) and 5G6 Detection (TA700126) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 ELISA with 3H5 Capture ([TA600125]) and 5G6 Detection (TA700126) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 Luminex with 2F9 Capture ([TA600123]) and 5G6 Detection (TA700126) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])



SOX17 Luminex with 3H5 Capture ([TA600125]) and 5G6 Detection (TA700126) Antibodies. Substrate used: Recombinant Human SOX17 ([TP320888])