

Product datasheet for TA700175

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ERG Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI8A9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8A9

Applications: ELISA, LMNX **Recommended Dilution:** 1:250 - 1:1000

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ERG(NP_891548) produced in HEK293T

cell.

Formulation: Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol

Concentration: 0.5 mg/ml **Conjugation:** Biotin

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Gene Name: ETS transcription factor ERG

Database Link: NP 891548

Entrez Gene 13876 MouseEntrez Gene 170909 RatEntrez Gene 2078 Human

P11308

Synonyms: erg-3; p55


Matched ELISA Pair: TA600174

Protein Families: Druggable Genome, Transcription Factors







Product images:



ERG Elisa with 5F12 Capture ([TA600174]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])




ERG Elisa with 4H7 Capture ([TA600175]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])




ERG Elisa with 5D10 Capture ([TA600178]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])






ERG Elisa with 4G10 Capture ([TA600179]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])




ERG Luminex Elisa with 5F12 Capture ([TA600174]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])




ERG Luminex Elisa with 4H7 Capture ([TA600175]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])





ERG Luminex Elisa with 5D10 Capture ([TA600178]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])



ERG Luminex Elisa with 4G10 Capture ([TA600179]) and 8A9 Detection (TA700175) Antibodies. Substrate used: Recombinant Human ERG ([TP308093])