

## Product datasheet for **TA700228**

### **TYRO3 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI4C8]**


#### **Product data:**

|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Clone Name:           | OTI4C8  |
| Applications:         | LMNX  |
| Recommended Dilution: | 1:250 - 1:1000  |
| Reactivity:           | Human, Mouse, Rat   |
| Host:                 | Mouse   |
| Isotype:              | IgG1  |
| Clonality:            | Monoclonal  |
| Immunogen:            | Full length human recombinant protein of human TYRO3 (NP_006284) produced in HEK293T cell.  |
| Formulation:          | Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol   |
| Concentration:        | 0.5 mg/ml   |
| Conjugation:          | Biotin  |
| Storage:              | Store at -20°C as received.   |
| Stability:            | Stable for 12 months from date of receipt.  |
| Gene Name:            | TYRO3 protein tyrosine kinase   |
| Database Link:        | <a href="#">NP_006284</a><br><a href="#">Entrez Gene 22174 Mouse</a> <a href="#">Entrez Gene 25232 Rat</a> <a href="#">Entrez Gene 7301 Human</a><br><a href="#">Q06418</a> |
| Synonyms:             | BYK; Dtk; Etk-2; Rek; RSE; Sky; Tif   |
| Matched ELISA Pair:   | TA600228  |
| Protein Families:     | Druggable Genome, Protein Kinase  |




[View online »](#)

**Product images:**



TYRO3 Luminex Elisa with 2G2 Capture ([TA600228]) and 4C8 Detection (TA700228) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])



TYRO3 Luminex Elisa with 9B12 Capture ([TA600229]) and 4C8 Detection (TA700228) Antibodies. Substrate used: Recombinant Human TYRO3 ([TP308260])