

Product datasheet for **TA700263**

Cystatin S (CST4) Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) **[Clone ID: OTI2C2]**


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2C2
Applications:	LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CST4 (NM_001899) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	cystatin S
Database Link:	NP_001890 Entrez Gene 1472 Human P01036
Synonyms:	MGC71923
Matched ELISA Pair:	TA600263



[View online »](#)

Product images:



CST4 Luminex Elisa with 2H10 Capture ([TA600263]) and 2C2 Detection (TA700263) Antibodies. Substrate used: Recombinant Human CST4 ([LY400707])