

Product datasheet for **TA700284**

CD19 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI2B11]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2B11
Applications:	LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD19 (NM_001770) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD19 molecule
Database Link:	NP_001761 Entrez Gene 930 Human P15391
Synonyms:	B4; CVID3
Matched ELISA Pair:	TA600286
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	B cell receptor signaling pathway, Hematopoietic cell lineage, Primary immunodeficiency



[View online »](#)

Product images:



CD19 Luminex Elisa with 2D3 Capture ([TA600286]) and 2B11 Detection (TA700284) Antibodies. Substrate used: Recombinant Human CD19 ([TP302922])