

## Product datasheet for **TA700313**

### **CD1C Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI2F4]**


#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI2F4
<b>Applications:</b>	ELISA, LMNX
<b>Recommended Dilution:</b>	1:250 - 1:1000
<b>Reactivity:</b>	Human
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human CD1C (NM_001765) produced in HEK293T cell.
<b>Formulation:</b>	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
<b>Concentration:</b>	0.5 mg/ml
<b>Conjugation:</b>	Biotin
<b>Storage:</b>	Store at -20°C as received.
<b>Stability:</b>	Stable for 12 months from date of receipt.
<b>Gene Name:</b>	CD1c molecule
<b>Database Link:</b>	<a href="#">NP_001756</a> <a href="#">Entrez Gene 911 Human P29017</a>
<b>Synonyms:</b>	BDCA1; CD1; CD1A; R7
<b>Matched ELISA Pair:</b>	TA600313, TA600314, TA600316
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	Hematopoietic cell lineage




[View online »](#)


**Product images:**




CD1C Elisa with 1F4 Capture ([TA600313]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])




CD1C Elisa with 2B3 Capture ([TA600314]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])




CD1C Elisa with 4C7 Capture ([TA600316]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 1F4 Capture ([TA600313]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 2B3 Capture ([TA600314]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 4C7 Capture ([TA600316]) and 2F4 Detection (TA700313) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])