

Product datasheet for **TA700318**

CD1C Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI2F1]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2F1
Applications:	ELISA, LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD1C (NM_001765) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD1c molecule
Database Link:	NP_001756 Entrez Gene 911 Human P29017
Synonyms:	BDCA1; CD1; CD1A; R7
Matched ELISA Pair:	TA600317
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage




[View online »](#)

Product images:



CD1C Elisa with 5A8 Capture ([TA600317]) and 2F1 Detection (TA700318) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])



CD1C Luminex Elisa with 5A8 Capture ([TA600317]) and 2F1 Detection (TA700318) Antibodies. Substrate used: Recombinant Human CD1C ([TP318490])