

Product datasheet for **TA700363**

POGK Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI3D2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI3D2
Applications:	ELISA
Recommended Dilution:	ELISA 1:250-1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human POGK (NM_017542) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	pogo transposable element derived with KRAB domain
Database Link:	NP_060012 Entrez Gene 57645 Human Q9P215
Synonyms:	BASS2; KRBOX2; LST003
Matched ELISA Pair:	TA600350
Protein Families:	Transcription Factors




[View online »](#)

Product images:



POGK Elisa with 2D5 Capture ([TA600350]) and 3D2 Detection (TA700363) Antibodies. Substrate used: Recombinant Human POGK ([TP322119])



POGK Luminex Elisa with 2D5 Capture ([TA600350]) and 3D2 Detection (TA700363) Antibodies. Substrate used: Recombinant Human POGK ([TP322119])