

Product datasheet for **TA700384**

p63 (TP63) Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI2G10]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2G10
Applications:	ELISA, LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human TP63 (NM_003722) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tumor protein p63
Database Link:	NP_003713 Entrez Gene 246334 Rat Entrez Gene 8626 Human Q9H3D4
Synonyms:	AIS; B(p51A); B(p51B); EEC3; KET; LMS; NBP; OFC8; p40; p51; p53CP; p63; p73H; p73L; RHS; SHFM4
Matched ELISA Pair:	TA600384
Protein Families:	Druggable Genome, Transcription Factors




[View online »](#)

Product images:



TP63 ELISA with 5A5 Capture ([TA600384]) and 2G10 Detection (TA700384) Antibodies. Substrate used: Recombinant Human TP63 ([TP308013])



TP63 Luminex Elisa with 5A5 Capture ([TA600384]) and 2G10 Detection (TA700384) Antibodies. Substrate used: Recombinant Human TP63 ([TP308013])