

Product datasheet for **TA700389**

Vimentin (VIM) Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI1A9]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1A9
Applications:	LMNX
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human VIM (NM_003380) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	vimentin
Database Link:	NP_003371 Entrez Gene 22352 MouseEntrez Gene 81818 RatEntrez Gene 7431 Human P08670
Synonyms:	CTRCT30; HEL113
Matched ELISA Pair:	TA600389
Protein Families:	ES Cell Differentiation/IPS



[View online »](#)

Product images:



VIM Luminex Elisa with 13B1 Capture ([TA600389]) and 1A9 Detection (TA700389) Antibodies. Substrate used: Recombinant Human VIM ([TP301546])