

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA700455

MAGEA3 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI24B6]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI24B6
Applications:	ELISA
Reactivity:	Human
Host:	Mouse
lsotype:	lgG2a
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human MAGEA3 (NP_005353) produced in HEK293T cell.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Predicted Protein Size:	34.6 kDa
Gene Name:	MAGE family member A3
Database Link:	<u>NP_005353</u> <u>Entrez Gene 4102 Human</u> <u>P43357</u>
Background:	Human Cancer/testis antigen MAGEA3 is a member of Melanoma Antigen Gene (MAGE) family, which has restricted expression to the testis and is aberrantly expressed in cancer cells. MAGEA3 has been found to be broadly expressed in a variety of malignancies such as melanoma, breast cancer, head and neck cancer, lung cancer, gastric cancer, skin squamous cell carcinoma, prostate cancer, colorectal cancer, etc. It has been shown that MAGEA3 expression is associated with cancer progression and metastasis, and increased expression of MAGE-A3 in cancer cells is associated with poor prognosis of cancer patients.
Synonyms:	CT1.3; HIP8; HYPD; MAGE3; MAGEA6



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US **MAGEA3** Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI24B6] – TA700455

Matched ELISA Pair: TA600453

Product images:



MAGEA3 ELISA with 1G9 Capture ([TA600453]) and 24B6 Detection (TA700455) Antibodies. Standard protein used: Recombinant Human MAGEA3 (HEK293T expressed protein Cat# [TP303288])

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US