

Product datasheet for **TA700458**

Galectin 3 (LGALS3) Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) **[Clone ID: OTI1A2]**


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI1A2
Applications:	ELISA
Recommended Dilution:	ELISA 1:250-1:1000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human LGALS3 (NM_002306) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	galectin 3
Database Link:	NP_002297 Entrez Gene 16854 Mouse Entrez Gene 3958 Human P17931
Synonyms:	CBP35; GAL3; GALBP; GALIG; L31; LGALS2; MAC2
Matched ELISA Pair:	TA600458
Protein Families:	Secreted Protein



[View online »](#)

Product images:



LGALS33 Luminex ELISA with 1A4 Capture ([TA600458]) and 1A2 Detection (TA700458) Antibodies. Substrate used: Recombinant Human LGALS33 ([TP308785])