

## Product datasheet for **TA700464**

### CD5 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OT15D4]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OT15D4
Applications:	ELISA
Recommended Dilution:	1:250 - 1:1000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human CD5 (NM_014207) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	CD5 molecule
Database Link:	<a href="#">NP_055022</a> <a href="#">Entrez Gene 921 Human</a> <a href="#">P06127</a>
Synonyms:	LEU1; T1
Matched ELISA Pair:	TA600464, TA600468
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage




[View online »](#)

Product images:



CD5 ELISA with 6F7 Capture ([TA600468]) and 5D4 Detection (TA700464) Antibodies. Substrate used: Recombinant Human CD5 ([TP306494])



CD5 ELISA with 10H3 Capture ([TA600464]) and 5D4 Detection (TA700464) Antibodies. Substrate used: Recombinant Human CD5 ([TP306494])