

Product datasheet for **TA700517**

PFDN3 (VBP1) Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) **[Clone ID: OTI2C2]**


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI2C2
Applications:	ELISA
Recommended Dilution:	ELISA 1:250-1:1000
Reactivity:	Human, Mouse
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human VBP1 (NM_003372) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	VHL binding protein 1
Database Link:	NP_003363 Entrez Gene 22327 Mouse Entrez Gene 7411 Human P61758
Synonyms:	HIBBJ46; PFD3; PFDN3; VBP-1
Matched ELISA Pair:	TA600050
Protein Families:	Druggable Genome




[View online »](#)

Product images:



VBP1 ELISA with Anti-Tag 1B12 Capture ([TA600050]) and 2C2 Detection (TA700517) Antibodies. Substrate used: Recombinant Human VBP1 ([TP308482])



VBP1 Luminex with Anti-Tag 1B12 Capture ([TA600050]) and 2C2 Detection (TA700517) Antibodies. Substrate used: Recombinant Human VBP1 ([TP308482])