

Product datasheet for **TA700525**

IKZF3 Biotinylated Mouse Monoclonal Detection Antibody (Biotin conjugated) [Clone ID: OTI5H2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI5H2
Applications:	ELISA
Recommended Dilution:	ELISA 1:250-1:1000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human IKZF3 (NM_012481) produced in HEK293T cell.
Formulation:	Stored in PBS (pH 7.4) with 0.05% sodium azide, 10mg/ml BSA, 50% glycerol
Concentration:	0.5 mg/ml
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	IKAROS family zinc finger 3
Database Link:	NP_036613 Entrez Gene 22780 Mouse Entrez Gene 303511 Rat Entrez Gene 22806 Human Q9UKT9
Synonyms:	AIO; AIOLOS; ZNFN1A3
Matched ELISA Pair:	TA600525




[View online »](#)

Product images:



IKZF3 ELISA with 2B9 Capture ([TA600525]) and 5H2 Detection (TA700525) Antibodies. Substrate used: Recombinant Human IKZF3 ([TP307547])



IKZF3 Luminex with 2B9 Capture ([TA600525]) and 5H2 Detection (TA700525) Antibodies. Substrate used: Recombinant Human IKZF3 ([TP307547])