

Product datasheet for **TA800775**


SOAT 2 (SOAT2) Mouse Monoclonal Antibody [Clone ID: OTI2D3]

Product data:


| | |
|-------------------------|--|
| Product Type: | Primary Antibodies |
| Clone Name: | OTI2D3 |
| Applications: | IHC |
| Recommended Dilution: | IHC 1:150 |
| Reactivity: | Human |
| Host: | Mouse |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Immunogen: | Human recombinant protein fragment corresponding to amino acids 1-120 of human SOAT2 (NP_003569) produced in E.coli. |
| Formulation: | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. |
| Concentration: | 1 mg/ml |
| Purification: | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Predicted Protein Size: | 59.7 kDa |
| Gene Name: | sterol O-acyltransferase 2 |
| Database Link: | NP_003569 Entrez Gene 8435 Human O75908 |
| Synonyms: | ACACT2; ACAT2; ARGP2 |
| Protein Families: | Transmembrane |
| Protein Pathways: | Steroid biosynthesis |




[View online »](#)

Product images:


Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human colon tissue using anti-SOAT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA800775)




Immunohistochemical staining of paraffin-embedded Human liver tissue within the normal limits using anti-SOAT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA800775)



Immunohistochemical staining of paraffin-embedded Carcinoma of Human lung tissue using anti-SOAT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA800775)



Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-SOAT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA800775)



Immunohistochemical staining of paraffin-embedded Human tonsil within the normal limits using anti-SOAT2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA800775)