

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA801095AM

# S100 beta (S100B) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2H7]

### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI2H7
Applications:	WB
Recommended Dilution:	WB 1:200
Reactivity:	Human, Mouse, Rat
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human S100B (NP_006263) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	10.5 kDa
Gene Name:	S100 calcium binding protein B
Database Link:	<u>NP_006263</u> Entrez Gene 20203 MouseEntrez Gene 25742 RatEntrez Gene 6285 Human P04271



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


# S100 beta (S100B) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI2H7] – TA801095AM

Background:The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-<br/>hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a<br/>wide range of cells, and involved in the regulation of a number of cellular processes such as<br/>cell cycle progression and differentiation. S100 genes include at least 13 members which are<br/>located as a cluster on chromosome 1q21; however, this gene is located at 21q22.3. This<br/>protein may function in Neurite extension, proliferation of melanoma cells, stimulation of<br/>Ca2+ fluxes, inhibition of PKC-mediated phosphorylation, astrocytosis and axonal<br/>proliferation, and inhibition of microtubule assembly. Chromosomal rearrangements and<br/>altered expression of this gene have been implicated in several neurological, neoplastic, and<br/>other types of diseases, including Alzheimer's disease, Down's syndrome, epilepsy,<br/>amyotrophic lateral sclerosis, melanoma, and type I diabetes. [provided by RefSeq, Jul 2008]

Synonyms:

NEF; S100; S100-B; S100beta

## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY S100B ([RC200277], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-S100B. Positive lysates [LY416735] (100ug) and [LC416735] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US