

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA801378AM

NM23A (NME1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI24A2]

Product data:

Product Type:	Primary Antibodies		
Clone Name:	OTI24A2		
Applications:	IHC, WB		
Recommended Dilution:	WB 1:2000, IHC 1:150		
Reactivity:	Human, Mouse, Rat		
Host:	Mouse		
lsotype:	lgG2b		
Clonality:	Monoclonal		
Immunogen:	Full length human recombinant protein of human NME1 (NP_937818) produced in E.coli.		
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.		
Concentration:	0.5 mg/ml		
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)		
Conjugation:	Biotin		
Storage:	Store at -20°C as received.		
Stability:	Stable for 12 months from date of receipt.		
Predicted Protein Size:	19.5 kDa		
Gene Name:	NME/NM23 nucleoside diphosphate kinase 1		
Database Link:	<u>NP_937818</u> Entrez Gene 18102 MouseEntrez Gene 191575 RatEntrez Gene 4830 Human P15531		




	NM23A (NME1) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI24A2] – TA801378AM	
Background:	This gene (NME1) was identified because of its reduced mRNA transcript levels in highly metastatic cells. Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by this gene) and 'B' (encoded by NME2) isoforms. Mutations in this gene have been identified in aggressive neuroblastomas. Two transcript variants encoding different isoforms have been found for this gene. Co-transcription of this gene and the neighboring downstream gene (NME2) generates naturally-occurring transcripts (NME1-NME2), which encodes a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Jul 2008]	
Synonyms:	AWD; GAAD; NB; NBS; NDKA; NDPK-A; NDPKA; NM23; NM23-H1	
Protein Families	: Druggable Genome, Stem cell - Pluripotency	
Protein Pathwa	Metabolic pathways, Purine metabolism, Pyrimidine metabolism	

wal Antikash. (Diatin saniumatash) (Clana


Product images:


170	_	
130	_	
100	_	
70	_	
55	_	
40	_	
35	_	
25	_	-
15	_	-
10	_	

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NME1 ([RC220517], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NME1. Positive lysates [LY404982] (100ug) and [LC404982] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])

Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])

Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])

Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])

Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-NME1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])





Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-NME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA801378])