

## Product datasheet for **TA801872AM**

### **MUC16 Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI3C6]**


#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI3C6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids around 13200-13700 of human MUC16 (NP_078966) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	mucin 16, cell surface associated
Database Link:	<a href="#">NP_078966</a> <a href="#">Entrez Gene 94025 Human</a> <a href="#">Q8WXI7</a>
Synonyms:	CA125




[View online »](#)

**Product images:**



E.coli lysate (left lane) and E.coli lysate expressing Human recombinant protein fragment (right lane) corresponding to amino acids 10329-10628 of human MUC16 (NP\_078966) were separated by SDS-PAGE and immunoblotted with anti-MUC16 monoclonal antibody (1:1000).



Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-MUC16 monoclonal antibody. 1:500`2h