

Product datasheet for TA802243BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OXSM Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1H6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1H6
Applications: IHC, WB

Reactivity: WB 1:2000, IHC 1:150 **Reactivity:** Human, Mouse, Rat

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 78-343 of human OXSM

(NP_060367) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 48.7 kDa

Gene Name: 3-oxoacyl-ACP synthase, mitochondrial

Database Link: NP 060367

Entrez Gene 71147 MouseEntrez Gene 289934 RatEntrez Gene 54995 Human

<u>Q9NWU1</u>

Background: This gene encodes a beta-ketoacyl synthetase. The encoded enzyme is required for

elongation of fatty acid chains in the mitochondria. Alternatively spliced transcript variants

have been described. [provided by RefSeq, Feb 2009]

Synonyms: FASN2D; KASI; KS


Protein Families: Druggable Genome






Protein Pathways: Fatty acid biosynthesis, Metabolic pathways


Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY OXSM (Cat# [RC201500], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OXSM(Cat# [TA802243]). Positive lysates [LY413479] (100ug) and [LC413479] (20ug) can be purchased separately from OriGene.





Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802243])



Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802243])







Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802243])

Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802243])