

Product datasheet for **TA802417S**

CD5 Mouse Monoclonal Antibody [Clone ID: OTI8E7]


Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Clone Name: | OTI8E7 |
| Applications: | FC |
| Recommended Dilution: | FLOW 1:50 |
| Reactivity: | Human |
| Host: | Mouse |
| Isotype: | IgG1 |
| Clonality: | Monoclonal |
| Immunogen: | Full length human recombinant protein of human CD5 (NP_055022) produced in HEK293T cell. |
| Formulation: | PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. |
| Concentration: | 1 mg/ml |
| Purification: | Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G) |
| Conjugation: | Unconjugated |
| Storage: | Store at -20°C as received. |
| Stability: | Stable for 12 months from date of receipt. |
| Predicted Protein Size: | 54.4 kDa |
| Gene Name: | CD5 molecule |
| Database Link: | NP_055022 Entrez Gene 921 Human P06127 |
| Synonyms: | LEU1; T1 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Hematopoietic cell lineage |




[View online »](#)


Product images:



Living HEK293T cells transfected with either [RC206494] overexpress plasmid (Red), compared to an IgG isotype control, (Green) or empty vector control plasmid (Blue) were immunostained by anti-CD5 antibody ([TA802417]), and then analyzed by flow cytometry (1:100).



Flow cytometric analysis of living 50ng/ml IL2-stimulated Jurkat cells, using anti-CD5 antibody ([TA802417]), (red), compared to an IgG isotype control, (green), and negative control (PBS), (blue) (1:100).



Flow cytometric analysis of living CCRF-CEM cells, using anti-CD5 antibody ([TA802417]), (red), compared to an IgG isotype control, (green), and negative control (PBS), (blue) (1:100).