

## Product datasheet for TA802432BM

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Ki67 (MKI67) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI13B11]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI13B11

Applications: IHC

Recommended Dilution: IHC 1:100~500

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1160-1493 of human

MKI67 (NP\_002408) produced in E.coli.

**Formulation:** PBS (pH 7.3) containing 1% BSA, 50% glycerol.

**Concentration:** 0.5 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

**Storage:** Store at -20°C as received.

**Stability:** Stable for 12 months from date of receipt.

**Gene Name:** marker of proliferation Ki-67

Database Link: NP 002408

Entrez Gene 4288 Human

P46013

**Background:** This gene encodes a nuclear protein that is associated with and may be necessary for cellular

proliferation. Alternatively spliced transcript variants have been described. A related

pseudogene exists on chromosome X. [provided by RefSeq, Mar 2009]


**Synonyms:** KIA; MIB-; MIB-1; PPP1R105

**Protein Families:** Druggable Genome, ES Cell Differentiation/IPS







## **Product images:**



Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-MKI67 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802432])




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-MKI67 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802432])



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-MKI67 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, [TA802432])





Immunohistochemical staining of paraffinembedded Mouse colon tissue using anti-MKI67 mouse monoclonal antibody. ([TA802432]; heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min)