

Product datasheet for TA803377BM

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GBA Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4G4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4G4
Applications: IHC, WB

Recommended Dilution: WB 1:500, IHC 1:150

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 40-315 of human GBA

(NP_000148) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol.

Concentration: 0.5 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: HRP

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 55.5 kDa

Gene Name: glucosylceramidase beta

Database Link: NP 000148

Entrez Gene 14466 MouseEntrez Gene 2629 Human

P04062

Background: This gene encodes a lysosomal membrane protein that cleaves the beta-glucosidic linkage of

glycosylceramide, an intermediate in glycolipid metabolism. Mutations in this gene cause

Gaucher disease, a lysosomal storage disease characterized by an accumulation of

glucocerebrosides. A related pseudogene is approximately 12 kb downstream of this gene on

chromosome 1. Alternative splicing results in multiple transcript variants. [provided by

RefSeq, Jan 2010]






Synonyms: GBA1; GCB; GLUC


Protein Families: Druggable Genome

Protein Pathways: Lysosome, Metabolic pathways, Other glycan degradation, Sphingolipid metabolism

Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY GBA (Cat# [RC216061], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GBA(Cat# [TA803377]). Positive lysates [LY400055] (100ug) and [LC400055] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 1 cell line by using anti-GBA monoclonal antibody (1:200).






Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA803377])




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA803377])



Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA803377])





Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA803377])