

Product datasheet for **TA803676BM**

IRAK2 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI1B4]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI1B4
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 362-625 of human IRAK2 (NP_001561) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	interleukin 1 receptor associated kinase 2
Database Link:	NP_001561 Entrez Gene 3656 Human O43187
Background:	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB. [provided by RefSeq, Jul 2008]
Synonyms:	IRAK-2
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: Apoptosis, Neurotrophin signaling pathway

Product images:



E.coli lysate (left lane) and E.coli lysate expressing human recombinant protein fragment corresponding to amino acids 362-625 of human IRAK2 (NP_001561) were separated by SDS-PAGE and immunoblotted with anti-IRAK2.