


## Product datasheet for **TA803760S**

### **NUDT10 Mouse Monoclonal Antibody [Clone ID: OTI1A6]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI1A6
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human NUDT10 (NP_694853) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	18.3 kDa
Gene Name:	nudix hydrolase 10
Database Link:	<a href="#">NP_694853</a> <a href="#">Entrez Gene 170685 Human</a> <a href="#">Q8NFP7</a>
Synonyms:	APS2; DIPPP3-alpha; DIPPP3a

[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NUDT10 ([RC208717], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NUDT10. Positive lysates [LY407127] (100ug) and [LC407127] (20ug) can be purchased separately from OriGene.