

Product datasheet for **TA803871AM**

MURF3 (TRIM54) Mouse Monoclonal Antibody (Biotin conjugated) [Clone ID: OTI7F8]


Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7F8
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 166-400 of human TRIM54 (NP_115935) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Biotin
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	tripartite motif containing 54
Database Link:	NP_115935 Entrez Gene 57159 Human Q9BYV2
Background:	The protein encoded by this gene contains a RING finger motif and is highly similar to the ring finger proteins RNF28/MURF1 and RNF29/MURF2. In vitro studies demonstrated that this protein, RNF28, and RNF29 form heterodimers, which may be important for the regulation of titin kinase and microtubule-dependent signal pathways in striated muscles. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]
Synonyms:	MURF; MURF-3; muRF3; RNF30




[View online »](#)


Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY TRIM54 (Cat# [RC218075], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TRIM54(Cat# [TA803871]). Positive lysates [LY403171] (100ug) and [LC403171] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffin-embedded Human Kidney tissue within the normal limits using anti-TRIM54 mouse monoclonal antibody. ([TA803871]) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Carcinoma of Human liver tissue using anti-TRIM54 mouse monoclonal antibody. (TA803871) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Carcinoma of Human lung tissue using anti-TRIM54 mouse monoclonal antibody. (TA803871) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Human pancreas tissue within the normal limits using anti-TRIM54 mouse monoclonal antibody. (TA803871) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Human endometrium tissue within the normal limits using anti-TRIM54 mouse monoclonal antibody. ([TA803871]) Dilution: 1:150




Immunohistochemical staining of paraffin-embedded Adenocarcinoma of Human endometrium tissue using anti-TRIM54 mouse monoclonal antibody. ([TA803871]) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Human prostate tissue within the normal limits using anti-TRIM54 mouse monoclonal antibody. ([TA803871]) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Carcinoma of Human prostate tissue using anti-TRIM54 mouse monoclonal antibody. (TA803871) Dilution: 1:150



Immunohistochemical staining of paraffin-embedded Human lymph node tissue within the normal limits using anti-TRIM54 mouse monoclonal antibody. (TA803871) Dilution: 1:150