

## Product datasheet for **TA804031BM**

### CD34 Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI51A1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI51A1
Applications:	FC, IF
Recommended Dilution:	FLOW 1:50
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	HEK293T cells transfected with full length human CD34 gene (RC204446)
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	40.5 kDa
Gene Name:	CD34 molecule
Database Link:	<a href="#">NP_001020280</a> <a href="#">Entrez Gene 947 Human</a> <a href="#">P28906</a>
Background:	The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Synonyms:	CD34 antigen; CD34 molecule; hematopoietic progenitor cell antigen CD34; OTTHUMP00000034733; OTTHUMP00000034734



[View online »](#)

**Protein Families:** Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transmembrane

**Protein Pathways:** Cell adhesion molecules (CAMs), Hematopoietic cell lineage

**Product images:**





Figure from citation: Immunofluorescence of CD34 protein level by using anti-CD34 antibody in human pulmonary arterial endothelial cells. Dilution: 1:100 [View Citation](#)



Living stable-transfected HEK293T cell line with [RC204446] overexpress plasmid (Red) or HEK293T cells (Blue) were immunostained by anti-CD34 antibody ([TA804031]), and then analyzed by flow cytometry.