

Product datasheet for **TA804192BM**


XPB (ERCC2) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI4B11]

Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI4B11
Applications:	WB
Recommended Dilution:	WB 1:2000
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 78-307 of human ERCC2 (NP_000391) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Gene Name:	ERCC excision repair 2, TFIIH core complex helicase subunit
Database Link:	NP_000391 Entrez Gene 2068 Human P18074
Synonyms:	COFS2; EM9; TFIIH; TTD; TTD1; XPB
Protein Families:	Druggable Genome, Transcription Factors
Protein Pathways:	Nucleotide excision repair



[View online »](#)

Product images:

Human recombinant protein fragment corresponding to amino acids 78-307 of human ERCC2 (NP_000391) produced in E.coli.