

Product datasheet for **TA804581**


FLIP (CFLAR) Mouse Monoclonal Antibody [Clone ID: OTI2H8]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI2H8
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:150
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 1-272 of human CFLAR (NP_003870) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	55.2 kDa
Gene Name:	CASP8 and FADD like apoptosis regulator
Database Link:	NP_003870 Entrez Gene 8837 Human O15519
Synonyms:	c-FLIP; c-FLIPL; c-FLIPR; c-FLIPS; CASH; CASP8AP1; Casper; CLARP; FLAME; FLAME-1; FLAME1; FLIP
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Apoptosis




[View online »](#)

Product images:


HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY CFLAR (Cat# [RC200608], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CFLAR antibody (Cat# TA804581). Positive lysates [LY418376] (100ug) and [LC418376] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (30ug per lane) from A549 lysate by using anti-FLIP monoclonal antibody (TA804581, 1:500).



Immunohistochemical staining of paraffin-embedded Carcinoma of Human prostate tissue using anti-CFLAR mouse monoclonal antibody. Heat-induced epitope retrieval by EDTA solution buffer pH 8.0 at 120°C for 3 min.