

## Product datasheet for **TA804638M**

### Angiopoietin 1 (ANGPT1) Mouse Monoclonal Antibody [Clone ID: OTI7F9]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7F9
Applications:	WB
Recommended Dilution:	WB 1:500
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 198-498 of human ANGPT1 (NP_001137) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57.3 kDa
Gene Name:	angiopoietin 1
Database Link:	<a href="#">NP_001137</a> <a href="#">Entrez Gene 11600 Mouse</a> <a href="#">Entrez Gene 89807 Rat</a> <a href="#">Entrez Gene 284 Human</a> <a href="#">Q15389</a>
Synonyms:	AGP1; AGPT; ANG1; HAE5
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)

## Product images:



Human recombinant protein fragment  
corresponding to amino acids 198-498 of human  
ANGPT1 (NP\_001137) produced in E.coli.