

## Product datasheet for **TA804679BM**

### Angiopoietin 1 (ANGPT1) Mouse Monoclonal Antibody (HRP conjugated) [Clone ID: OTI8H2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI8H2
Applications:	IHC, WB
Recommended Dilution:	WB 1:2000, IHC 1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 198-498 of human ANGPT1 (NP_001137) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol.
Concentration:	0.5 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	HRP
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	57.3 kDa
Gene Name:	angiopoietin 1
Database Link:	<a href="#">NP_001137</a> <a href="#">Entrez Gene 11600 Mouse</a> <a href="#">Entrez Gene 89807 Rat</a> <a href="#">Entrez Gene 284 Human</a> <a href="#">Q15389</a>



[View online »](#)

**Background:**


Angiopoietins are proteins with important roles in vascular development and angiogenesis. All angiopoietins bind with similar affinity to an endothelial cell-specific tyrosine-protein kinase receptor. The protein encoded by this gene is a secreted glycoprotein that activates the receptor by inducing its tyrosine phosphorylation. It plays a critical role in mediating reciprocal interactions between the endothelium and surrounding matrix and mesenchyme and inhibits endothelial permeability. The protein also contributes to blood vessel maturation and stability, and may be involved in early development of the heart. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Dec 2010]

**Synonyms:**


AGP1; AGPT; ANG1; HAE5

**Protein Families:**

Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein

**Product images:**

Human recombinant protein fragment corresponding to amino acids 198-498 of human ANGPT1 (NP\_001137) produced in E.coli.



Immunohistochemical staining of paraffin-embedded Carcinoma of Human liver tissue using anti-ANGPT1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA804679])