

## Product datasheet for **TA804793**


### Calpain 5 (CAPN5) Mouse Monoclonal Antibody [Clone ID: OTI7F2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OTI7F2
Applications:	IHC
Recommended Dilution:	IHC 1:150
Reactivity:	Human, Mouse, Rat
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 274-578 of human CAPN5 (NP_004046) produced in E.coli.
Formulation:	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Conjugation:	Unconjugated
Storage:	Store at -20°C as received.
Stability:	Stable for 12 months from date of receipt.
Predicted Protein Size:	73 kDa
Gene Name:	calpain 5
Database Link:	<a href="#">NP_004046</a> <a href="#">Entrez Gene 12337 Mouse</a> <a href="#">Entrez Gene 171495 Rat</a> <a href="#">Entrez Gene 726 Human</a> <a href="#">O15484</a>
Synonyms:	ADNIV; HTRA3; nCL-3; VRNI
Protein Families:	Druggable Genome, Protease



[View online »](#)

**Product images:**

Immunohistochemical staining of paraffin-embedded Carcinoma of Human liver tissue using anti-CAPN5 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804793)