

Product datasheet for TA804826S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CLOCK Mouse Monoclonal Antibody [Clone ID: OTI2G10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2G10
Applications: IHC, WB

Recommended Dilution: WB 1:2000, IHC 1:150

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-285 of human CLOCK

(NP_004889) produced in E.coli.

Formulation: PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Conjugation: Unconjugated

Storage: Store at -20°C as received.

Stability: Stable for 12 months from date of receipt.

Predicted Protein Size: 95.1 kDa

Gene Name: clock circadian regulator

Database Link: NP 004889

Entrez Gene 12753 MouseEntrez Gene 60447 RatEntrez Gene 9575 Human

<u>O15516</u>

Synonyms: bHLHe8; KAT13D


Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: Circadian rhythm - mammal






Product images:



Human recombinant protein fragment corresponding to amino acids 1-285 of human CLOCK (NP_004889) produced in E.coli.



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-CLOCK mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, [TA804826])